

## CORNER TABLES - :L Shape, U Shape:-

METROLABS Laboratory Tables are modular type. The modular will be designed in such a way that the modules can be dismantled very easily in case of maintenance if required.

The under storage modules will be independent and entire load on the work surface will be sustained by the "C" frame / Pedestal frame / Base frame and not on the storage modules. Under bench storage modules will have drawer and bottom cupboard and front open shutter.

The drawer and cupboards will have necessary locking arrangements. Components to be powder coated is pre-treated using seven tank processor, to ensure that the substrate is clean and free from grease, rust. Powder coating thickness is 80 micron.

